Millipore®

Preparation, Separation,
Filtration & Monitoring Products

MERCK

FLASH® Allergen Indicator Protein Test


Cleaning Verification for Protein Removal and Allergen Control Programs.

FLASH® is a total protein visual test that rapidly detects protein residues left on food contact surfaces after cleaning. Protein is a difficult food residue to remove. Most allergens are proteins, so quick verification of surface hygiene helps minimize the risk of cross-contact to allergen-free products. Using FLASH® regularly also helps reduce the opportunity for biological contamination.

As part of a HACCP allergen management program, FLASH® supports process verification requirements that ensure cleaning methods, independently validated to effectively remove allergens, are consistently applied.

Used together with the MVP ICON® for ATP testing, both can be incorporated into validated cleaning procedures (SSOPs) as part of a more comprehensive HACCP management program.

To find out more, visit:
SigmaAldrich.com

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada

A perfect complement to ATP testing


Two Applications. You Decide.

FLASH® is a self-contained sampling device containing a pre-moistened tip and all necessary reagents. One swab provides two sensitivity levels, allowing customers to choose the appropriate level based on application- or area-specific limits of detection for protein residue.

Manage Risks. Be Prepared. Everytime.

- Detects total protein including “Big 8” food allergens
- Rapid pass / fail result within 10 minutes
- Single product, multi-sensitivity detection. (20 µg at room temp / 3 µg when incubated)
- Simple online training available 24/7 with certification for audit compliance

3 - Step Process. Simple.

SWAB

Sample a 4 x 4 inch
(10 x 10 cm) surface

ACTIVATE

Activate swab and
shake for 5 seconds

READ


Use chart to
interpret result
(Pass, Caution or Fail)

To place an order or receive technical assistance in the U.S. and Canada, call toll-free 1-800-245-0113
For other countries across Europe and the world, please visit: SigmaAldrich.com

© 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

Merck, the vibrant M, FLASH, and MVP ICON are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

2018-FLASH-R000_AS_BR_MK MK_BR3081EN

FLASH[®] RT Allergen Indicator Total Protein Test

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Validation and Application of FLASH Rapid Protein Detection Test

FLASH® is a rapid total protein detection test which serves as verification for adequate removal of proteinaceous matter, including allergens, on manufacturing surfaces.

Monitoring adequate removal after cleaning allows for immediate corrective action and reduces the potential risk for cross-contamination of product and food safety incidents, especially during change-overs.

FLASH is a self-contained sampling device containing a pre-moistened tip and all necessary reagents. It is designed with two levels of sensitivity to accommodate varying industry needs. When used at room temperature, FLASH can detect as little as 20 µg of total protein within 10 minutes. If greater sensitivity is needed, such as prior to processing critical products or when sampling

key manufacturing areas, FLASH can be incubated for 10 ± 4 minutes in a dry heat bath at 70 °C for detection down to 3 µg of total protein.

Sensitivity Validation & Performance

In an in-house study, FLASH was tested against common allergenic proteins covering the Food Safety and Inspection Service's (FSIS) "Big 8", as well as against total proteins such as raw ground beef and raw ground turkey. Refer to Table 1 for study results, which indicate the approximate level where tests began showing a positive (Fail) result. All proteins tested were detected by FLASH with an average detection level of 20 µg at room temperature and 3 µg at 70 °C.

Results Interpretation

FLASH is a qualitative visual test interpreted via color change. All tests are initially green to start. The solution will remain green if no protein is detected. The solution will turn an increasing shade of purple with increasing amounts of protein detected.

Table 1 – FLASH Detection of Proteins in Foods

Protein
Milk Powder
Raw Ground Turkey
Soy Flour
Peanut Butter
Egg Powder
Raw Fish (Cod)
Roasted Almonds
Raw Shrimp
Raw Ground Beef
Sesame
Gluten Flour

If the solution turns grey, trace amounts of protein are detectable. The surface should be re-cleaned and then re-tested per SSOP's before proceeding.

Simple To Use:

SWAB

Sample a 4 x 4 in (10 x 10 cm) surface

ACTIVATE

Activate swab and shake for 5 seconds

READ

Use chart to interpret result (Pass, Caution or Fail)

BIOCONTROL
www.biocontrolsys.com

Results. Right Now.

North America: tel +1.425.603.1123 bcs_us@biocontrolsys.com
 Asia Pacific: tel +852 5808.2443 bcs_asiapacific@biocontrolsys.com
 Brazil: tel +55 19.3213.5919 bcs_latnamerica@biocontrolsys.com
 France: tel +33 (0)4.7271.5680 bcs_fr@biocontrolsys.com
 Italy: tel +39 06.9148.831 bcs_it@biocontrolsys.com
 Netherlands: tel +31 (0)180.333.955 bcs_nl@biocontrolsys.com
 United Kingdom: tel +44 (0)845.539.9902 bcs_uk@biocontrolsys.com