

Biotech

&

Scientific

THINK OF LAB .THINK OF US  
คิดถึงแลป คิดถึงเรา

---

**CATALOG**  
**BIOTECHSCI**®

---

**Biotech and Scientific.co.th**  
**155 onnut 17 yek 16 Suanglaung**  
**Suanglaung Bangkok 10250 Thailand**  
**Tel +66 2 7177915 Fax +66 2 717 7919**  
**E.mail: biotechsci@hotmail.com**  
**Mobile : 0859171259**



Biotechsci®

1000mL  
BIO32006

1000  
800  
600  
400  
200  
BORO 3  
APPROX VOL

BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



### Bottles, Reagent Clear Screw Neck

- Autoclavable PP screw cap and pouring ring
- These complete bottles are autoclavable & can be sterilized
- These bottles are chemical resistance
- Printed with bio code.
- Clear scale Graduation & large marking area




BiotechSciCODE	Item Code	Capacity (ml)	Diameter (mm)	Height (mm)	Pack of	
BIO32001		25*	36	70	10	GL32
BIO32002	01-VK15-01	50*	45	86.5	10	GL32
BIO32003	01-VK15-02	100	56	100	10	GL45
BIO32004	01-VK15-03	250	69.5	138	10	GL45
BIO32005	01-VK15-04	500	90	175.5	10	GL45
BIO32006	01-VK15-05	1000	104.5	225	10	GL45
BIO32007	01-VK15-06	2000	134.5	260	10	GL45
BIO32008		5000	181	330	1	GL45
BIO32009		10000	227	410	1	GL45
BIO32010		20000	288	505	1	GL45

### Bottles, Reagent Amber Screw Neck

- Autoclavable PP screw cap and pouring ring
- These complete bottles are autoclavable and can be sterilized
- These bottles are chemical resistance
- Printed with bio code.
- Clear scale Graduation & large marking area



BiotechSci CODE	Item Code	Capacity (ml)	sub (ml)	Diameter (mm)	Height (mm)	Pack of	
BIO32001A		25*		36	70	10	GL32
BIO32002A		50*		46	87	10	GL32
BIO32003A	01-VK46-02	100	10	60.5	114	10	GL45
BIO32004A		250		70	138	10	GL45
BIO32005A	01-VK46-04	500	50	88.5	175.5	10	GL45
BIO32006A	01-VK46-05	1000	100	105	227	10	GL45
BIO32007A		2000		136	260	10	GL45
BIO32008A		5000		181	330	1	GL45
BIO32009A		10000		227	410	1	GL45
BIO32010A		20000		288	505	1	GL45





### Beakers, Low form with Spout

- Designed with pouring spout
- Clear Graduation for easy reading convenience
- wall thickness
- Large labeling field for easy marking
- Ideal for heating liquids
- Highly resistant to chemical



BiotechSci CODE	ITEMCODE	Capacity(ml)	Height (mm)	Diameter (mm)	Packof
BIO13001		5	30	22	10
BIO13002		10	35	26	10
BIO13003		25	50	34	10
BIO13004	01-VK14-01	50	57.5	41	10
BIO13005	01-VK14-10	100	70	50	10
BIO13006	01-VK14-12	150	88	60	10
BIO13007	01-VK14-02	250	94.5	69	10
BIO13008	01-VK14-03	400	115.5	80	10
BIO13009	01-VK14-11	500	115	83.5	8
BIO13010	01-VK14-04	600	126.5	90.5	8
BIO13011		800	135	100	6
BIO13012	01-VK14-05	1000	157	105	6
BIO13013	01-VK14-06	2000	188.5	131	4
BIO13014	01-VK14-07	3000	210	150	1
BIO13015	01-VK14-08	5000	260	172.5	1
BIO13016	01-VK14-09	10000	350	220	1

## Beakers, Tall form with Spout

- Designed with pouring spout
- Clear Graduation for easy reading convenience
- wall thickness
- Large labeling field for easy marking
- Ideal for heating liquids
- Highly resistant to chemical




BiotechSciCODE	Item Code	Capacity(ml)	Dia (mm)	Height(mm)	Pack of
BIO17001	01-VK16-01	50	38	70	10
BIO17002	01-VK16-02	100	48	80	10
BIO17003		150	54	95	10
BIO17004	01-VK16-03	250	60	120	10
BIO17005	01-VK16-04	400	70	132	10
BIO17006	01-VK16-05	500	75	142	8
BIO17007		600	80	150	8
BIO17008	01-VK16-06	1000	90.5	190.5	6
BIO17009	01-VK16-07	25	34	64.5	10
BIO17010	01-VK16-08	300	64	118	12






## Flask, Erlenmeyer with Joint

- Flask ,Conical with Joint
- Body is thickwall



BiotechSciCODE	ItemCode	capacity(ml)	SocketSize	H.(mm)	D.(mm)	Pack of
BIO38001		25	14/23	70	42	10
BIO38002		25	19/26	70	42	10
BIO38003		50	14/23	85	51	10
BIO38004		50	19/26	85	51	10
BIO38005		50	24/29	85	51	10
BIO38006		50	29/32	85	51	10
BIO38007		100	14/23	100	64	10
BIO38008		100	19/26	100	64	10
BIO38009		100	24/29	100	64	10
BIO38010		100	29/32	100	64	10
BIO38011		150	19/26	118	72	10
BIO38012		150	24/29	118	72	10
BIO38013		150	29/32	118	72	10
BIO38014		200	29/32	131	79	10
BIO38015		250	19/26	85	10	10
BIO38016	01-VK40-01	250	24/29	145	85	10
BIO38017	01-VK40-02	250	29/32	143	85.5	10
BIO38018		300	29/32	156	87	10
BIO38019		500	19/26	175	105	10
BIO38020	01-VK40-04	500	24/29	174	105	10
BIO38021	01-VK40-03	500	29/32	175	105	10
BIO38022	01-VK40-05	1000	24/29	220	130.5	10
BIO38023	01-VK40-06	1000	29/32	221	128	10
BIO38024		1000	34/35	220	131	10
BIO38025		2000	29/32	270	166	6
BIO38026		2000	34/35	270	166	6
BIO38027		3000	45/50	310	187	1
BIO38028		5000	45/50	365	220	1





BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259


E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



**Flask, Round Bottom, Narrow Neck,**

**- Manufactured borosilicate glass**

**- wall thickness**



BiotechSci CODE	ItemCODE	capacity(ml)	D.(mm)	D1.(mm)	H.(mm)	Packof
BIO43001		50	51	26	90	10
BIO43002	01-VK24-09	100	64	26	110	10
BIO43003	01-VK24-01	250	84.5	34	141.5	10
BIO43004		500	105	34	168	10
BIO43005		1000	131	42	200	6
BIO43006	01-VK24-05	2000	166	50	240	6
BIO43007	01-VK24-06	3000	185	50	260	1
BIO43008	01-VK24-07	5000	223	50	305	1
BIO43009	01-VK24-08	10000	279	65	380	1
BIO43010		20000	350	75	510	1
BIO43011	01-VK24-03	125	69	28	128.5	10
BIO43012	01-VK24-04	150	77	30	130	10

BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




**Flat Bottom, Narrow Mouth**

- borosilicate glass 3.3

- wall thickness




BiotechSci CODE	Item Code	capacity (ml)	D. (mm)	D1. (mm)	H. (mm)	Pack of
BIO45001		50	51	26	90	10
BIO45002		100	64	26	110	10
BIO45003		250	85	34	143	10
BIO45004	01-VK25-01	500	105	34	168	10
BIO45005		1000	131	42	200	6
BIO45006	01-VK25-03	2000	166	50	240	6
BIO45007	01-VK25-04	3000	185	50	260	1
BIO45008	01-VK25-05	5000	223	50	305	1
BIO45009	01-VK25-06	10000	279	65	380	1
BIO45010		20000	350	75	510	1
BIO45011	01-VK25-02	125	68	28.5	120.5	10



**Round Bottom Flask, Single Neck with Joint**

BiotechSci	Item Code	capacity	socket	Height	Diameter	Pack
BIO23000		25 ML	14/23	85	41	10
BIO23001		25 ML	19/26	85	41	10
BIO23002		25 ML	24/29	85	41	10
BIO23003		50 ML	14/23	90	51	10
BIO23004		50 ML	19/26	90	51	10
BIO23005		50 ML	24/29	90	51	10
BIO23006		50 ML	29/32	90	51	10
BIO23038		100 ML	14/23	105	64	10
BIO23039		100 ML	19/26	105	64	10
BIO23007	01-VK02-12	100 ML	24/29	105	64	10
BIO23008	01-VK02-13	100 ML	29/32	105	64	10
BIO23009		150 ML	19/26	137	75	10
BIO23010		150 ML	24/29	137	75	10
BIO23040		150 ML	29/32	137	75	10
BIO23041		250 ML	14/23	140	85	10
BIO23042		250 ML	19/26	140	85	10
BIO23011	02-VK02-17	250 ML	24/29	140	85	10
BIO23012	02-VK02-18	250 ML	29/32	140	85	10
BIO23013		250 ML	34/35	140	85	10
BIO23014		500 ML	19/26	163	105	10
BIO23015	01-VK02-22	500 ML	24/29	163	105	10
BIO23016	01-VK02-23	500 ML	29/32	163	105	10
BIO23017		500 ML	34/35	163	105	10
BIO23018		1000 ML	19/26	200	131	10
BIO23019	01-VK02-06	1000 ML	24/29	200	131	10
BIO23020	01-VK02-07	1000 ML	29/32	200	131	10
BIO23021		1000 ML	34/35	200	131	10
BIO23022		2000 ML	24/29	240	165	6
BIO23023		2000 ML	29/32	240	165	6
BIO23024		2000 ML	34/35	240	165	6
BIO23025		3000 ML	24/29	260	185	1
BIO23026		3000 ML	29/32	260	185	1
BIO23027		3000 ML	34/35	260	185	1
BIO23028		3000 ML	45/50	260	185	1
BIO23029		5000 ML	24/29	305	223	1
BIO23030		5000 ML	29/32	305	223	1
BIO23031		5000 ML	34/35	305	223	1
BIO23032		5000 ML	45/40	305	223	1
BIO23033		10000 ML	29/32	380	279	1
BIO23034		10000 ML	34/35	380	279	1
BIO23035		10000 ML	45/40	380	279	1
BIO23036		10000 ML	55/44	380	279	1
BIO23037		20000 ML	55/44	435	345	1



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259


E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)

---



### **Volumetric flasks Class “A”**

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- Calibrated to contain (In) Class “A” standards.
- Volumetric flasks are supplied with polypropylene stopper.




## Volumetric flasks Class "A"

BiotechSci CODE	Item CODE	Volume	Tolerance	Height	Diameter	Neck size	Quantity in box
		*1 ml	±0.025 ml	65 mm	7.50 mm	NS 07/16	2 pieces
		*2 ml	±0.025 ml	70 mm	7.50 mm	NS 07/16	2 pieces
		*3 ml	±0.025 ml	70 mm	7.50 mm	NS 07/16	2 pieces
		*4 ml	±0.025 ml	70 mm	7.50 mm	NS 07/16	2 pieces
		5 ml	±0.025 ml	70 mm	7.50 mm	NS 07/16	2 pieces
BIO1001	01-VK11-05	5 ml	±0.025 ml	70 mm	10.00 mm	NS 10/19	2 pieces
BIO1002		10 ml	±0.025 ml	90 mm	10.00 mm	NS 10/19	2 pieces
BIO1003		20 ml	±0.040 ml	110 mm	10.00 mm	NS 10/19	2 pieces
BIO1004	01-VK11-06	25 ml	±0.030 ml	110 mm	10.00 mm	NS 10/19	2 pieces
BIO1004A	01-VK10-02	25 ml	±0.040 ml	110 mm	12.50 mm	NS 12/21	3 pieces
BIO1005		25 ml	±0.040 ml	110 mm	12.50 mm	NS 12/20	2 pieces
BIO1005A		25 ml					
BIO1006	01-VK11-01	50 ml	±0.060 ml	140 mm	12.50 mm	NS 12/21	2 pieces
BIO1007		50 ml	±0.060 ml	140 mm	14.50 mm	NS 14/23	2 pieces
BIO1008		100 ml	±0.100 ml	170 mm	12.50 mm	NS 12/21	2 pieces
BIO1009	01-VK11-02	100 ml	±0.100 ml	170.50 mm	14.50 mm	NS 14/23	2 pieces
BIO1009A	01-VK10-01	100 ml	±0.100 ml	170.50 mm	14.50 mm	NS 14/23	2 pieces
BIO1010		*150 ml	±0.150 ml	190 mm	14.50 mm	NS 14/21	2 pieces
BIO1011		200 ml	±0.150 ml	210 mm	14.50 mm	NS 14/22	2 pieces
BIO1012	01-VK11-03	250 ml	±0.150 ml	220 mm	14.50 mm	NS 14/23	2 pieces
BIO1013		*300 ml	±0.200 ml	230 mm	14.50 mm	NS 14/23	2 pieces
BIO1014		*400 ml	±0.200 ml	250 mm	18.80 mm	NS 19/26	2 pieces
BIO1015	01-VK11-04	500 ml	±0.250 ml	260 mm	19.00 mm	NS 19/26	2 pieces
BIO1016		1000 ml	±0.400 ml	300 mm	24.00 mm	NS 24/29	2 pieces
BIO1017		2000 ml	±0.600 ml	370 mm	29.20 mm	NS 29/32	2 pieces
BIO1018		5000 ml	±1.200 ml	475 mm	34.50 mm	NS 34/35	1 piece
BIO1019		*10000 ml	±2.000 ml	570 mm	45.00 mm	NS 45/40	1 piece

## Burettes

*Clear burettes* PTFE stopcock

- Manufactured from calibrated borosilicate 3.3 glass tubes.
- Calibrated to the tolerances are according to DIN/ISO Class "AS" standards.



Item CODE	volume	tolerance	Subdivision	quantity in box
-----------	--------	-----------	-------------	-----------------

BIO4001	10 ml	±0.020 ml	±0.020 ml	1 piece
BIO4001A	10 ml	±0.020 ml	±0.020 ml	1 piece
BIO4002	01-VK47-01 25 ml	±0.030ml	±0.030 ml	1 piece
BIO4002A	01-VK47-02 25 ml	±0.030ml	±0.030 ml	1 piece
BIO4003	01-VK47-03 50 ml	±0.050 ml	±0.050 ml	1 piece
BIO4003A	01-VK47-04 50 ml	±0.050 ml	±0.050 ml	1 piece

*A is Scale Amber colour*





## Automatic burettes class AS

### *Clear & amber automatic burettes class "AS"*

- Manufactured from borosilicate glass.
- VALVE stopcock and also an intermediate PTFE stopcock

BiotechSci CODE	Item CODE	volume	tolerance	subdivision	packing	in box
BIO5001		10 ml	±0.020 ml	0.02 ml	1 piece	
BIO5001A		10 ml	±0.020 ml	0.02 ml	1 piece	
BIO5002	01-VK22-01	25 ml	±0.030 ml	0.05 ml	1 piece	
BIO5002A	01-VK27-01	25 ml	±0.030 ml	0.05 ml	1 piece	
BIO5003	01-VK22-02	50 ml	±0.050 ml	0.10 ml	1 piece	
BIO5003A	01-VK27-02	50 ml	±0.050 ml	0.10 ml	1 piece	

\*\*\*A is Amber



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



---

### ***Clear & amber schilling burettes class "AS"***

**- Described as automatic burettes with fast titration and practical usage.**

**- 10 ml burette has a polypropylene reservoir bottle of 500 ml volume capacity where as 25 and 50 ml burettes are used with bottles of 1000 ml volume capacity.**

<i>BiotechSci CODE</i>	<i>Item CODE</i>	<i>volume</i>	<i>tolerance</i>	<i>subdivision</i>
<i>BIO6001</i>		<i>10 ml</i>	<i>±0.030 ml</i>	<i>0.02 ml</i>
<i>BIO6001A</i>	<i>01-VK44.01</i>	<i>10 ml</i>	<i>±0.030 ml</i>	<i>0.02 ml</i>
<i>BIO6002</i>	<i>01-VK43.01</i>	<i>25 ml</i>	<i>±0.050 ml</i>	<i>0.10 ml</i>
<i>BIO6002A</i>	<i>01-VK44.02</i>	<i>25 ml</i>	<i>±0.050 ml</i>	<i>0.10 ml</i>
<i>BIO6003</i>	<i>01-VK43.02</i>	<i>50 ml</i>	<i>±0.050 ml</i>	<i>0.10 ml</i>
<i>BIO6003A</i>	<i>01-VK44.03</i>	<i>50 ml</i>	<i>±0.050 ml</i>	<i>0.10 ml</i>

**High form measuring cylinders class "A"**

- Manufactured from borosilicate glass 3.3  
which is resistant to heat and almost all chemicals.  
Calibrated to contain (In) Class "A" standards.

**Thick and durable plastic base of the cylinder**


**does not have the risk of breakage**

**during application The plastic base fits the cylinder**

**body perfectly well and can be easily removed.**



BiotechSci CODE	Item Code	Volume	Tolerance	Subdivision	Height	Diameter	Quantity in box
BIO2001	01-VK12-02	5 ml	±0.050 ml	±0.100 ml	115 mm	14 mm	2 pieces
BIO2002		10 ml	±0.100ml	±0.200ml	140 mm		2 pieces
BIO2003		25 ml	±0.250 ml	±0.500 ml	170 mm		2 pieces
BIO2004		50 ml	±0.500 ml	±1.000 ml	200 mm		2 pieces
BIO2005	01-VK12-01	100 ml	±0.500 ml	±1.000 ml	260 mm	32.00 mm	2 pieces
BIO2006	01-VK12-03	250 ml	±1.000 ml	±2.000 ml	335 mm	42.5 mm	2 pieces
BIO2007		500 ml	±2.000 ml	±5.000 ml	390 mm		2 pieces
BIO2008		1000 ml	±5.000 ml	±10.000 ml	470 mm		2 pieces
BIO2009		2000 ml	±10.000 ml	±20.000 ml	570 mm		1 piece



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



## *Mixing cylinders class "A"*

Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.

- Calibrated to contain (In) Class "A" standards.
- polypropylene stopper.

BiotechSci CODE	Item CODE	volume	Tolerance	Height	neck size	quantity in box
BIO3001		10 ml	±0.100 ml	160 mm	NS 10/19	2 pieces
BIO3002		25 ml	±0.250ml	195 mm	NS 14/23	2 pieces
BIO3003		50 ml	±0.500 ml	230 mm	NS 19/26	2 pieces
BIO3004	01-VK31-01	100 ml	±0.500 ml	290 mm	NS 24/29	2 pieces
BIO3005		250 ml	±1.000 ml	370 mm	NS 29/32	2 pieces
BIO3006		500 ml	±2.500 ml	395 mm	NS 34/35	2 pieces
BIO3007		1000 ml	±5.000 ml	510 mm	NS 45/40	1 piece



### Filter flasks glass - glass side arm

- Manufactured from borosilicate glass 3.3 .
- Filter flasks are manufactured with glass side arms
- 6-8 mm hose diameter.




BiotechSci CODE	Item Code	volume	Diameter top	Diameter neck	Diameter bottom	Hieght	quantity in box
BIO29001	01-VK05-01	100 ml		24 mm			1 piece
BIO29002	01-VK05-03	250 ml	40 mm	28 mm	85.5 mm	153 mm	1 piece
BIO29003	01-VK05-04	500 ml	34.5 mm	34 mm		186 mm	1 piece
BIO29004	01-VK05-05	1000 ml		45 mm	135 mm	225 mm	1 piece
BIO29005	01-VK05-06	2000 ml	51 mm	60 mm	162 mm	300 mm	1 piece
BIO29006	01-VK05-07	5000 ml		100 mm			1 piece

### Filtering set glass 47 mm complete set

- Manufactured from borosilicate glass 3.3 .
- With Glass Connector
- Made from heavy wall which can withstand high vacuum
- Printed in white color enamel without Graduation
- Filter Diameter, mm: 47
- Funnel Capacity, mL: Funnel: 300 ml Flask:1 L
- Outlet Fitting: 6 mm (1/4 in.) O.D. tubulated cap sidearm to vacuum



BiotechSci CODE	ItemCode	Detail	Pack
BIO31001	01-VK33-01	filtering set 47mm	1 set



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




## TEST TUBE WITHOUT CAP BOROSILIGATE GLASS

BiotechSci CODE	Item Code	O.D.xHeight(mm)	Pack of
BIO12001	01-vk18-02	16x125mm	100
BIO12002	01-vk18-01	16x100mm	100
BIO12003	01-vk18-03	16x150mm	100
BIO12004	01-vk18-04	18x150mm	100
BIO12005	01-vk18-05	25x150mm	100
BIO12006	01-vk18-06	25x200mm	100
BIO12007	01-vk18-07	20x125mm	100
BIO12008	01-vk18-08	18x125mm	100
BIO12009	01-vk18-09	22x100mm	100
BIO12010	01-vk18-10	22x125mm	100
BIO12011	01-vk18-11	22x150mm	100
BIO12012	01-vk18-12	22x200mm	100
BIO12013	01-vk18-13	25x175mm	100
BIO12014	01-vk18-14	32x200mm	100
BIO12015	01-vk18-15	38x200mm	100



## TEST TUBE WITH SCREW CAP BOROSILIGATE GLASS

BiotechSci CODE	Item Code	O.D. x Height (mm)	Pack of
BIO11012		32x200mm	100
BIO11013	01-vk19-13	38x200mm	100
BIO11014		13x100mm	100
BIO11001	01-vk19-01	16x125mm	100
BIO11002	01-vk19-02	16x150mm	100
BIO11003	01-vk19-03	20x150mm	100
BIO11004	01-vk19-04	20x125mm	100
BIO11005	01-vk19-05	25x150mm	100
BIO11006	01-vk19-06	22x175mm	100
BIO11007	01-vk19-07	18x150mm	100
BIO11008	01-vk19-08	25x175mm	100
BIO11009	01-vk19-09	22x125mm	100
BIO11010	01-vk19-10	22x150mm	100
BIO11011	01-vk19-11	22x200mm	100
BIO11016	01-vk19-16	25x200mm	100
BIO11015	01-vk19-15	20x200mm	100



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



### Centrifuge tube without screw cap




BiotechSci CODE	Item Code	Capacity (ml)	Dia (mm)	Hight (mm)	pack of
BIO49001	01-VK39-01	15	15	121	20
BIO49002	01-VK39-02	50	25	180	10

### Centrifuge tube with screw cap



BiotechSci CODE	Item Code	Capacity (ml)	Dia (mm)	Hight 1 (mm)	pack of
BIO48001	01-VK41-02	15	15	121	20
BIO48002	01-VK41-01	50	25	180	10





BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



## separating funnels squibb

- Manufactured from borosilicate glass 3.3 which is resistant to heat and chemicals.
- PTFE stopcock.
- With SCALE



BiotechSci CODE	Item CODE	volume	Diameter	neck NS size	packsize	
BIO25001 w/o	BIO25001	50 ml		NS 19/26	1 piece	
BIO25002 w/o	BIO25002	100 ml		NS 19/26	1 piece	
BIO25003 w/o	BIO25003	01-VK35-01	250 ml	85 mm	NS 29/32	1 piece
BIO25004 w/o	BIO25004	01-VK35-02	500 ml	102 mm	NS 29/32	1 piece
BIO25005 w/o	BIO25005	01-VK35-03	1000 ml		NS 29/32	1 piece
BIO25006 w/o	BIO25006	01-VK35-04	2000 ml		NS 29/32	1 piece

BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



## Funnels glass


*Laboratory funnels glass*

*-Manufactured from*

*borosilicate glass 3.3*



BiotechSci	Item code	funnel diameter	stem diameter	stem length
BIO27001		40 mm	6 mm	40 mm
BIO27002	01-VK17-03	70 mm	7 mm	70 mm
BIO27003	01-VK17-01	80 mm	8 mm	80 mm
BIO27004		100 mm	10 mm	100 mm
BIO27005	01-VK17-04	120 mm	16 mm	120 mm
BIO27006	01-VK17-05	150 mm	16 mm	150 mm



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




## BUCHNER FUNNEL GLASS

\*\*\*\*\*BOROSIRIGATE GLASS\*\*\*\*\*

\*\*\*\*\*USE FOR FILLTRATION\*\*\*\*\*



BiotechSciCODE	ItemCODE	Capacity	ofDise	Porosity	Packof
BIO36001	01-VK34-01	35 ML	30	1	2
BIO36002	01-VK34-06	35 ML	30	2	2
BIO36003	01-VK34-07	35 ML	30	3	2
BIO36004	01-VK34-08	35 ML	30	4	2
BIO36005	01-VK34-02	80 ML	40	1	2
BIO36006	01-VK34-09	80 ML	40	2	2
BIO36007	01-VK34-10	80 ML	40	3	2
BIO36008	01-VK34-11	80 ML	40	4	2
BIO36009	01-VK34-03	200 ML	65	1	2
BIO36010	01-VK34-12	200 ML	65	2	2
BIO36011	01-VK34-13	200 ML	65	3	2
BIO36012	01-VK34-14	200 ML	65	4	2
BIO36013	01-VK34-04	500 ML	90	1	2
BIO36014	01-VK34-15	500 ML	90	2	2
BIO36015	01-VK34-16	500 ML	90	3	2
BIO36016	01-VK34-17	500 ML	90	4	2
BIO36017	01-VK34-05	1000 ML	120	1	2
BIO36018	01-VK34-18	1000 ML	120	2	2
BIO36019	01-VK34-19	1000 ML	120	3	2
BIO36020	01-VK34-20	1000 ML	120	4	2



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




## CRUSIBLE GLASS

- BOROSIRIGATE GLASS



BiotechSci CODE	Item Code	Porosity	Capacity	Dia of Dise mm	Height	Pack of
BIO35001	01-VK32-01	1	15 ML	20	50	10
BIO35002	01-VK32-04	2	15 ML	20	50	10
BIO35003	01-VK32-05	3	15 ML	20	50	10
BIO35004	01-VK32-06	4	15 ML	20	50	10
BIO35005	01-VK32-02	1	30 ML	30	60	10
BIO35006	01-VK32-07	2	30 ML	30	60	10
BIO35007	01-VK32-08	3	30 ML	30	60	10
BIO35008	01-VK32-09	4	30 ML	30	60	10
BIO35009	01-VK32-03	1	50 ML	40	65	10
BIO35010	01-VK32-10	2	50 ML	40	65	10
BIO35011	01-VK32-11	3	50 ML	40	65	10
BIO35012	01-VK32-12	4	50 ML	40	65	10
BIO35013	01-VK32-13	5	50 ML	40	65	10



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




## Desiccators glass & plastic

Superior quality desiccators  
are manufactured from  
borosilicate glass 3.3



BiotechSci CODE	Item CODE	desiccator diameter	quantity in box
BIO47001 nonvac		200 mm	1 piece
BIO47002 nonvac	01-VK20-02	250 mm	1 piece
BIO47003 nonvac	01-VK20-01	300 mm	1 piece
BIO47001 vac		200 mm	1 piece
BIO47002 vac		250 mm	1 piece
BIO47003 vac		300 mm	1 piece



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



*Cal*



*Cal+ Thermo*

## Density bottles

BiotechSciCODE	Item CODE	volume	tolurance
BIO9001 Uncal		5 ml	
BIO9001 Cal		5 ml	
BIO9001 Cal+Thermo		5 ml	
BIO9002 Uncal		10 ml	
BIO9002 Cal	01-VK37-01	10 ml	± 0.3 ml
BIO9002 Cal+Thermo		10 ml	
BIO9003 Uncal		25 ml	
BIO9003 Cal	01-VK37-02	25 ml	
BIO9003 Cal+Thermo	01-VK38-01	25 ml	± 0.8 ml
BIO9004 Uncal		50 ml	
BIO9004 Cal	01-VK37-03	50 ml	± 1.0 ml
BIO9004 Cal+Thermo	01-VK38-02	50 ml	± 1.0 ml

BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



## Dishes, Crystallizing, without Spout

- Crystallizing Dishes
- Made of Boro 3.3
- Have BIG LABEL
- Thinkness wall



BiotechSciCODE	Item Code	Ø (mm) x h	Capacity (ml)	Pack of
BIO33001	01-VK42-01	40x25 mm	60 ml	10
BIO33002	01-VK42-02	50x30 mm	65 ml	10
BIO33004	01-VK42-03	60x35 mm	80 ml	10
BIO33005	01-VK42-04	70x40 mm	180 ml	10
BIO33006	01-VK42-05	80x45 mm	200 ml	10
BIO33007	01-VK42-06	95x55 mm	280 ml	10
BIO33008	01-VK42-07	100x50 mm	300 ml	10
BIO33009	01-VK42-08	150x75 mm	1200 ml	10
BIO33010	01-VK42-09	190x100 mm	2500 ml	10
BIO33011	01-VK42-10	170x90 mm	1770 ml	10

BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259


E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



## Alihn Condenser

### -Borosiligate glass

BiotechSci	Item Code	SIZE	JOINT	Stopcock
BIO15001		150 MM	24/40	GLASS
BIO15002		160 MM	14/23	GLASS
BIO15003		160 MM	19/26	GLASS
BIO15004		250 MM	14/23	GLASS
BIO15005		250 MM	19/23	GLASS
BIO15006		250 MM	24/29	GLASS
BIO15007		250 MM	24/40	GLASS
BIO15008		300 MM	24/32	GLASS
BIO15009		300 MM	24/29	GLASS
BIO15010	01-VK30-01	300 MM	29/32	GLASS
BIO15011		400 MM	24/40	GLASS
BIO15012		400 MM	29/32	GLASS
BIO15013		400 MM	24/40	GLASS
BIO15014		500 MM	24/40	GLASS
BIO15015		150 MM	24/40	GL 14
BIO15016		160 MM	14/23	GL 14
BIO15017		160 MM	19/26	GL 14
BIO15018		250 MM	14/23	GL 14
BIO15019		250 MM	19/26	GL 14
BIO15020		250 MM	24/29	GL 14
BIO15021		250 MM	24/40	GL 14
BIO15022		300 MM	24/29	GL 14
BIO15023		300 MM	24/40	GL 14
BIO15024		300 MM	29/32	GL 14
BIO15025		400 MM	24/29	GL 14
BIO15026		400 MM	29/32	GL 14
BIO15027		400 MM	24/40	GL 14
BIO15028		500 MM	24/40	GL 14





BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259



E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



## Coli Condenser

### -BOROSILIGATE 3.3

BiotechSci	Item Code	SIZE	JOINT	Stopcock
BIO16001		150 MM	24/40	GLASS Connector
BIO16002		160 MM	14/23	GLASS Connector
BIO16003		160 MM	19/26	GLASS Connector
BIO16004		250 MM	14/23	GLASS Connector
BIO16005		250 MM	19/23	GLASS Connector
BIO16006		250 MM	24/29	GLASS Connector
BIO16007		250 MM	24/40	GLASS Connector
BIO16008		300 MM	24/32	GLASS Connector
BIO16009		300 MM	24/29	GLASS Connector
BIO16010	01-VK29-01	300 MM	29/32	GLASS Connector
BIO16011		400 MM	24/40	GLASS Connector
BIO16012		400 MM	29/32	GLASS Connector
BIO16013		400 MM	24/40	GLASS Connector
BIO16014		500 MM	24/40	GLASS Connector
BIO16015		150 MM	24/40	GL 14
BIO16016		160 MM	14/23	GL 14
BIO16017		160 MM	19/26	GL 14
BIO16018		250 MM	14/23	GL 14
BIO16019		250 MM	19/26	GL 14
BIO16020		250 MM	24/29	GL 14
BIO16021		250 MM	24/40	GL 14
BIO16022		300 MM	24/29	GL 14
BIO16023		300 MM	24/40	GL 14
BIO16024		300 MM	29/32	GL 14
BIO16025		400 MM	24/29	GL 14
BIO16026		400 MM	29/32	GL 14
BIO16027		400 MM	24/40	GL 14
BIO16028		500 MM	24/40	GL 14



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




## Soxhlet extractors

-BOROSILIGATE GLASS

HIGH HEAT RESISTANT



BiotechSci CODE	Item CODE	volume	cone NS neck	socket NS neck	quantity in box	
BIO24001 w/o		30 ml	NS 29/32	NS 29/32	1 piece	
BIO24002 w/o	BIO24002	70 ml	NS 29/32	NS 34/35	1 piece	
BIO24003 w/o	BIO24003	100 ml	NS 29/32	NS 45/40	1 piece	
BIO24004 w/o	BIO24004	150 ml	NS 29/32	NS 45/40	1 piece	
BIO24005 w/o	BIO24005	250 ml	NS 29/32	NS 45/40	1 piece	
BIO24006 w/o	BIO24006	01-VK10-04	500 ml	NS 29/32	NS 45/40	1 piece
BIO24007 w/o	BIO24007	01-VK10-05	1000 ml	NS 29/33	NS 45/41	1 piece
BIO24008 w/o	BIO24008	01-VK10-06	2000 ml	NS 29/34	NS 45/42	1 piece



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."


Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



ADAPTER, Receiver

BiotechSci CODE	Item Code	Socket	Approx Length mm	Pack
BIO21001	01-VK07-01	14/23	65	10
BIO21002	01-VK07-02	19/26	65	10
BIO21003	01-VK07-03	24/29	65.5	10
BIO21004	01-VK07-04	29/32	65.5	10
BIO21005	01-VK07-05	24/35	65	10
BIO21006	01-VK07-06	40/48	65	10
BIO21007	01-VK07-07	45/40	65.5	10
BIO21008	01-VK07-08	50/42	65	10
BIO21009	01-VK07-09	55/44	65	10



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




ADAPTER, Multiple 1 pararell and 1 joint

BiotechSci CODE	Item Code	Socket	Cone	Pack
BIO19001		14/23	14/23	10
BIO19002		14/23	19/26	10
BIO19003		19/26	19/26	10
BIO19004		19/26	24/29	10
BIO19005	01-VK06-04	24/29	24/29	10
BIO19006		19/26	34/35	10
BIO19007	01-VK06-06	29/32	29/32	10



ADAPTER, Bend for socket

BiotechSciCODE	ItemCode	Socket	Cone	Pack
BIO51000		14/23	14/23	1
BIO51001		24/29	14/23	1
BIO51002		19/26	19/26	1
BIO51003		24/29	19/26	1
BIO51004		29/32	19/26	1
BIO51006	01-VK08-04	29/32	29/32	1
BIO51005	01-VK08-03	24/29	24/29	1
BIO51007		34/35	24/29	1
BIO51008		29/32	24/29	1
BIO51009		34/35	34/35	1



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)




---

## Adapter Three way



BiotechSciCode	Item Code	Socket	Cone	Cone size to fit condensor	Vertical	Horizontal	Pack
BIO52001		14/23	14/23	14/23			1
BIO52002		14/23	19/26	19/26			1
BIO52003	01-VK09-05	24/29	24/29	24/29	159.5	129.5	1
BIO52004	01-VK09-06	29/32	29/32	29/32	147	103.5	1
BIO52005		14/23	34/35	19/26			1
BIO52006	01-VK09-03	14/23	24/29	24/29	144.5	95.5	1
BIO52007		14/23	34/35	24/29			1
BIO52008		14/23	29/32	24/29			1
BIO52009	01-VK09-04	14/23	29/32	29/32	155	124	1
BIO52010		14/23	34/35	29/32			1



BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB . THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : [Biotechsci@hotmail.com](mailto:Biotechsci@hotmail.com) ; [www.biotechsci.co.th](http://www.biotechsci.co.th)



## Watch glass

*-Made from*

*BOROSILIGATE 3.3*



BiotechSci CODE	Item CODE	Diameter	Thickness	pcs
BIO53001	01-VK45-01	50 mm	1.8 mm	12
BIO53002	01-VK45-02	60 mm	1.8 mm	12
BIO53003	01-VK45-03	75 mm	1.8 mm	12
BIO53004	01-VK45-04	100 mm	2.0 mm	12
BIO53005	01-VK45-05	125 mm	2.2 mm	12
BIO53006	01-VK45-06	150 mm	2.2 mm	12