

Biotechsci®

Pipette

Volume Micropipette

121°C
FULLY
AUTOCLAVE

Easy Click Stop Volume Setting.

A pipettes which combine a fully autoclavable setup with a lightweight and ergonomic body for a hassle and pain-free pipetting experience while maintaining unmatched accuracy and precision..

Biotechsci code	Item code	Volume Range [µl] ปริมาตร	R ≤ [µl] % ค่าความถูกต้อง	CV ≤ [µl] % ค่าความถูกต้อง	Division [µl]
BIO54001	08-VK03-01	0.5 - 10	1.2 0.12	0.8 0.08	0.1
BIO54002	08-VK03-02	5 - 50	1.0 0.5	0.4 0.2	0.5
BIO54003	08-VK03-03	20 - 200	0.8 1.6	0.3 0.6	1.0
BIO54004	08-VK03-04	100 - 1000	0.8 8.0	0.3 3.0	5.0
BIO54005	08-VK03-05	1000 - 5000	0.8 40.0	0.3 15.0	50.0
BIO54006	08-VK03-06	1000 - 10000	0.8 50.0	0.3 15.0	50.0



Two Step Plunger Operation

Easy to read display

Fully Autoclavable.
Highest Resistance To Physical
And Chemical Corrosion

- Ergonomic and Light Weight Design.
- Smooth Plunger Mechanism Offers Least Resistance For Ease-of-use And Long Hours Of Pipetting.
- Durable Tip Cone Made Of PVDF (Polyvinylidene Fluoride) Offers The Highest Resistance To Physical And Chemical Corrosion And Is Compatible With Wide Range Of Tips.
- Tip Ejector For Smooth One-hand Ejection Of The Tip Preventing Any Contamination.
- Two Step Plunger Operation Allows Forward & Reverse Pipetting Technique.
- Individual calibration certificate according to ISO8655
- Easy calibration and Fully Autoclave
- Easy maintenance

BIOTECH AND SCIENTIFIC CO.LTD

BIOTECHSCI GLASSWARE LAB " THINK OF LAB THING US."

Tel +66 2 7177915 MoBile +66 8 5 9171259

E.mail : Biotechsci@hotmail.com ; www.biotechsci.co.th